

C17orf58-chromosome 17 open reading frame 58 Gene

Reference: PTXBC036760

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UPF0450 protein C17orf58; chromosome 17 open reading frame 58		
For research use only.			

GENE DESCRIPTION	chromosome 17 open reading frame 58
GENE ID	284018
ORIGIN SPECIES	Human
ACCESSIONS	BC036760
VECTOR	pDONR223
NCBI SYMBOL	C17orf58
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaccgctgtacctacccccgacggcttcttctccgagtcacatgttagccctggactcctccagctgcaataagccgtgtccagagtttaaacctggcagcaggtatattgtgatggccacatctaccataagagacggcagctccctacagctctgctccaggtcctgagagcgcgctccgtccaggggatggactgctaaggagcagcagcagctatgtgaaaagggttaatcgaaaagggaaggcagattcaaggtgcagttcatacccaatgcatttga