

ZNF396-zinc finger protein 396 Gene

Reference: PTXBC036765

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ZSCAN14; zinc finger protein 396; zinc finger and SCAN domain-containing protein 14		

For research use only.

GENE DESCRIPTION	zinc finger protein 396
GENE ID	252884
ORIGIN SPECIES	Human
ACCESSIONS	BC036765
VECTOR	pDONR223
NCBI SYMBOL	ZNF396
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctgcaaatgggaaagtcacatcactcctaacacaaactcagaggagtgtatgggattctgacagagaagatggagaggaagagcagacctgtgatccagactctagcctccactggagcagcagctacagccagagacctccgccgaattcaggcagtttggtacc
aggattcactgggccccatgaggctctgagccggtctggaacttftcatctctggtgagccggaagtgcacaccaaggagcagatcctggagctgctggtgctggagcagttcctggccatcctcccaaaagagctcaggcctgggtgcagaagcatccagagaatgga
gaggaaactgtgactatgctggaggatgtggagagagagcttgatggaccaagcaggttaagaaggatgcctgtggagatgaacccccaggctgagtaagagcactggcacatcattgggagtag