

CTF1-cardiotrophin 1 Gene

Reference: PTXBC036787

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CT-1; CT1; cardiotrophin-1; cardiophin 1		
For research use only.			

GENE DESCRIPTION	cardiotrophin 1
GENE ID	1489
ORIGIN SPECIES	Human
ACCESSIONS	BC036787
VECTOR	pDONR223
NCBI SYMBOL	CTF1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagccggaggagggaagctctggacccccagactgattctcagttcacttctccccacttggaggccaagatccgtcagacacacagccttgcgcacctctccaccaataacgtgagcagctgctccaggaaatgtgagctccaggagaccccttcgggtgccagcttct
cgccgcccggctgcccgtggccggcctgagcggcccggtccgagccacgcccgggtgccagtgacagagcggctgcccgtggagcggcgccgtgcccgcgtgccccgctgctggagcagtggtgcccagggccgagctgaacccgcccgcgcgcgcctgc
tgcgccgctggaggacgcccgcgcagcccccggccctggggccgcccgtggaggccttctgcccgcgtggcgccgccaaccgcccggccggagacccccccgcccaccgctcagccgctccgccaccgggtcttccccgccaaggtgctggggctcc
gcgtttgcccctctaccgcagtggtgagccgaccgaggggcagctggccagctgctgcccggggctcggcctga