



Reference: PTXBC036789

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CPX2; UNQ676; inactive carboxypeptidase-like protein X2; 4632435C11Rik; carboxypeptidase Hlo; carboxypeptidase X (M14 family), member 2; carboxypeptidase-like protein X2; cytosolic carboxypeptidase; carboxypeptidase X, M14 family member 2		

GENE DESCRIPTION	carboxypeptidase X (M14 family), member 2
GENE ID	119587
ORIGIN SPECIES	Human
ACCESSIONS	BC036789
VECTOR	pDONR223
NCBI SYMBOL	CPXM2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]