

C12orf68-chromosome 12 open reading frame 68 Gene

Reference: PTXBC036801

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C12orf68; coiled-coil domain-containing protein 184; coiled-coil domain containing 184		

For research use only.

GENE DESCRIPTION	chromosome 12 open reading frame 68
GENE ID	387856
ORIGIN SPECIES	Human
ACCESSIONS	BC036801
VECTOR	pDONR223
NCBI SYMBOL	C12orf68
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggacggtctgctggagatcatgaccaaggacggcgacatgccggccccctggaggtgtccaccgtgccggcagtggtggggagtgatctccggggagtacaacggcgccatgaaggaaactgatggagcacctgaaagcccagctgcaagccctgtttgaggacgtga
gggcatgagggggccctggacgagcaggcctcgcacatccaggtgctctcggacgagtggtgcgccaaccagcagccatcgtctccatgtgccagattatgaccactgcgccccgccaggcgctgtggcggtgctggcgggcaaggggagcttcagagcgaccccaag
agccggagactccttcgctgggacgggacagcgcttctggtgctgcatcccgaggacgaggagcaggaagaagagaaggagatgccagccccgccacaccctcagtcactgtgagcggccgaaagccctgtgctgtctcttggggggacgggcccactgt
ggagcccctcgatatgccgacattaccctgctgcaactggaggcgagggcctcctgtga