

TMEM131-transmembrane protein 131 Gene

Reference: PTXBC036841

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CC28; PRO1048; RW1; YR-23; transmembrane protein 131; 2610524E03Rik		

For research use only.

GENE DESCRIPTION	transmembrane protein 131
GENE ID	23505
ORIGIN SPECIES	Human
ACCESSIONS	BC036841
VECTOR	pDONR223
NCBI SYMBOL	TMEM131
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagttctcttcttcttcttcttcaagtaccctaaataaaacagccaaatggaagcaaacacaagttgacaaggcagcctcgtccgggcaagaacggcaaccccaactttgtgcagtcacggctggctacgacaagagccaggtgggaatggccttgctaaagttcttcaaca
aaacaggtttctccagcagccttggcatttcacagctcgtgtgacagcgatggctcagacagctcgggtttgtggagtcctcagcaacccaagcagcctgacttcactccctcaattcgttctccgcttggaaactctttaaactaactgggtgaagtttcagcaaacctcgattatct
cgatcgtgcaatcaggcctcacagaggagctggaacgagtttaaatagtgcccttcataccttgggagtcgccaagcagatccagtccttctgtgccagcagttccgctcccgaccacacagccacatcggtctcgttaacaccagcggcctgtgtggtccaccactccatca
gcagctccatttgggtccagcaaccttagcagcgccttcccttaccactccagcaaacacagctggcaagcatggcctcatgggcacagaaaactccctgctcctcagctccctccacctccagtcagctgacgacttgggacagacctacaacccgtggcgatattggagcccca
cgattggaagaagaagctcgacccttggcttaattgcacttctcagcagaattaa