

C6orf141-chromosome 6 open reading frame 141 Gene

Reference: PTXBC036917

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C6orf141; chromosome 6 open reading frame 141		
For research use only.			

GENE DESCRIPTION	chromosome 6 open reading frame 141
GENE ID	135398
ORIGIN SPECIES	Human
ACCESSIONS	BC036917
VECTOR	pDONR223
NCBI SYMBOL	C6orf141
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagaccggggcctcagggagctgcgaatccatggactctccgcagcctgggggacctcgggcttttccgaggagtagggcggggctccgtagctccggcgccgggaatcccgacggcagggcgagccgaagccagggcgccgacaggac
agaacggcagatcgggcccctcggacctcggccggggaggaattggaccgtgagtcctgggtcagagagaaagtctcttctcctgcacccagagaggtggttagggactcagggggaccctgcacgggaagaggtggccggtgcagaggaccttctcatcggtggagagg
accacggcagggagggcccaactaccttctgtcttcaacgagaaaaagcgaattctggcaggcgtgtagccccgcgggacgcagacccacccaaatagctgttgcgggtggaggattatcaggtaacacaagaagtgttcagacctctggccaaagggtcgcagac
cacgaggactgaggagcacttctgacccgcgtcacttttcgagcagtagagagggacagcctggggaacgtggggccctgctgaatcccggtctccaggcagcaacaggggcatcccgctccacgcccggggaggggtttcacctctccagggacttggctcgag
gaaattaaactctaa