



Reference: PTXBC036923

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C9orf150; HYST0841; LRAP35b; bA3L8.2; leucine rich adaptor protein 1-like; leucine repeat adaptor protein 35b; leucine rich adaptor protein 1 like		

For research use only.

GENE DESCRIPTION	chromosome 9 open reading frame 150
GENE ID	286343
ORIGIN SPECIES	Human
ACCESSIONS	BC036923
VECTOR	pDONR223
NCBI SYMBOL	C9orf150
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]