



Reference: PTXBC036951

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CCDC52; SPICE; spindle and centriole-associated protein 1; coiled-coil domain-containing protein 52; spindle and centriole associated protein 1		

For research use only.

GENE DESCRIPTION	coiled-coil domain containing 52
GENE ID	152185
ORIGIN SPECIES	Human
ACCESSIONS	BC036951
VECTOR	pDONR223
NCBI SYMBOL	CCDC52
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]