

FLJ20433-hypothetical protein FLJ20433 Gene

Reference: PTXBC037355

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	mut-7; exonuclease mut-7 homolog; exonuclease 3'-5' domain-containing protein 3; testicular tissue protein Li 125; exonuclease 3'-5' domain containing 3		

For research use only.

GENE DESCRIPTION	hypothetical protein FLJ20433
GENE ID	54932
ORIGIN SPECIES	Human
ACCESSIONS	BC037355
VECTOR	pDONR223
NCBI SYMBOL	FLJ20433
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaccagagatcccgtgtgacccctgccgtggcgcagcgcaccgcacatggccaggacccctcctctgctcctgcaggccctgcagaccctgtgtccacgcgggagcggagcagctccggaggagcctggcgggggtttgctgcttggacgacccctggccggg
cttctggacatgctggagagctgccggggccaggcgggagaggccctccctggcggcctggatctccaccagctgcagtgctggctacaggcacagccatgccgagcctggccagcacagcctgaggctgaagcagctgcaggcccgagcgggtcaagtctcactgaga
gccccccagccttggcggcaccactggccagcactctcagctgcaggatgcagacaggagctgcctgtgacacagctccaccgcctccaccagaggcaggttcagagaagcagccacgctgggcacgacgttgaagctgcagtcggagcttggcgttgaagaatgagcaltc
ccactgctctccaggacaaggctggcctcgtggagcgtatgtggccgcttcccgacctccagaggagctgctgtcctcatggattcctgtgtgccagcccggctttgacatcaaggacgttgccagctgccatgtgactgcagctggagcagctgcctctggagaggagcggct
gggctgagcagagcacagccctcttctcccgagtgga