



Reference: PTXBC037567

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AGXT2L2; PHLU; 5-phosphohydroxy-L-lysine phospho-lyase; 5-phosphonooxy-L-lysine phospho-lyase; alanine--glyoxylate aminotransferase 2-like 2		

For research use only.

GENE DESCRIPTION	alanine-glyoxylate aminotransferase 2-like 2
GENE ID	85007
ORIGIN SPECIES	Human
ACCESSIONS	BC037567
VECTOR	pDONR223
NCBI SYMBOL	AGXT2L2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1