

RPS19BP1-ribosomal protein S19 binding protein 1 Gene

Reference: PTXBC037573

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AROS; S19BP; active regulator of SIRT1; 40S ribosomal protein S19-binding protein 1; RPS19-binding protein 1; homolog of mouse S19 binding protein; ribosomal protein S19 binding protein 1		

For research use only.

GENE DESCRIPTION	ribosomal protein S19 binding protein 1
GENE ID	91582
ORIGIN SPECIES	Human
ACCESSIONS	BC037573
VECTOR	pDONR223
NCBI SYMBOL	RPS19BP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccgccccctgtcgccggggcctggagctgctggcggcgccgagggccccgggacccctcaggtcaggccaagccgagaggggctccggtgaaacggccccggaagacgaaggcaattcaggcccagaactcgggaactcggccaagggaagggtgcccaagt
cggcactggagctaccggaagcgagagtgctgagaccctcagagtaaacctgaagtcttgaccaggacgagaagcaccgtggctgagctgtgagccagcagattttgcccagaaccggggccgcaaggcctgtgaccggcctgtgccaagaccaagaagaaggc
tgagggcaccgtgttcaccgaggaagacttcagaagttccagcagggaatacttcggcagctag