

YAF2-YY1 associated factor 2 Gene

Reference: PTXBC037777

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	YY1-associated factor 2; YY1 associated factor 2		

For research use only.

GENE DESCRIPTION	YY1 associated factor 2
GENE ID	10138
ORIGIN SPECIES	Human
ACCESSIONS	BC037777
VECTOR	pDONR223
NCBI SYMBOL	YAF2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

```
atgggagacaagaagagccccaccaggccgaagcggcagcccgaagccgtcctcggatgagggttactgggactgtagcgtctgcacctccggaacagcgcgagggccttcaagtgcattgatgtgcgatgtcggaaagggcacctccaccgatccacattgtttgaagtattgtga
gtgcctcagaggaaggaaccgttgaattccaattcaggaaggaacctgcacctgtctccagttgggtgcacagcaggttactcagcagttgtgcctctacacagtcacaagaaagagaaaaagataaagttagaaaaagaaaaagtgaaaaggaaacaactagcaaaaaaga
atagccataagaaaaccaggccaagattgaaaaatgtggatcggagtagtctcagcattggaaattactgttggagatctgacagtcattattacagacttaaggagaaaacaagtcaccgcctgcatctagtctgcctctgcagatcaacacagtcacaagcggctctagctctgataa
cacagagaggaatgtccaggtcatctcaccagaggagaagcctcatcattgaatggagaatctcattaa
```