



Reference: PTXBC037900

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C-terminal-binding protein 2; ribeye; C-terminal binding protein 2		
For research use only.			

GENE DESCRIPTION	C-terminal binding protein 2
GENE ID	1488
ORIGIN SPECIES	Human
ACCESSIONS	BC037900
VECTOR	pDONR223
NCBI SYMBOL	CTBP2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]