



Reference: PTXBC037906

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C4orf19; chromosome 4 open reading frame 19		
For research use only.			

GENE DESCRIPTION	chromosome 4 open reading frame 19
GENE ID	55286
ORIGIN SPECIES	Human
ACCESSIONS	BC037906
VECTOR	pDONR223
NCBI SYMBOL	C4orf19
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

atggggctgcagggtgctgtaaaaatacaaaagctatctcttggatcagtgcaagtgccctctcctggctatgcaatgaagtgcaacagctgcaagctagatgaagacgacactgataaaataaaaggcaaatggagcagtggaagctctggtcagaaaaatgacctcagaggcagggtcctaaagaagactgagacgacgacgacgactgagggcagctgctcctcctcagcaacaggccgctccacagggggagcctggagggagcaaacacatgctgctgtgtaacggctgagggccctctgctcactcacaacgaccttggaaattccagccctccagccagggatgacagggtctctggccagctatgccaactatgcttcccacaactcaacctctctggaaaggggccgacagaaacagactgctgctctgctgcctcagaaaggcccaactgaggaatggagactccagactcctctggagcagaaatttctctggagacaaacactgcttccagataccagccccagattacctcagcattggggccagctggagcaacgcttgatcataatgaaaaggactgtgtttcagaaccatactgaggatgaatccctgaggggaattcagccccacttgggggagcattggttgaatacgcctctctgtgaggaagctgggatctaatgaatgagatgtgaaacagaagttctaagcatctgtttaatgagaagggtctgttcatgagcctgtggttgactcaggaacagcagggagacacacatcggctcgtgagatgggaatggggagattgtgacgaggatgcagcgggtggcgaggcccttgcactttagaagctgctactgcagcagaagattggatgagcagattgaa