

ADAT2-adenosine deaminase, tRNA-specific 2, TAD2 homolog (S. cerevisiae) Gene

Reference: PTXBC037955

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	tRNA-specific adenosine-34 deaminase subunit ADAT2; DEADC1; TAD2; dJ20N2; dJ20N2.1; tRNA-specific adenosine deaminase 2; adenosine deaminase, tRNA-specific 2, TAD2 homolog; deaminase domain containing 1; deaminase domain-containing protein 1; tRNA-specifi		

For research use only.

GENE DESCRIPTION	adenosine deaminase, tRNA-specific 2, TAD2 homolog (S. cerevisiae)
GENE ID	134637
ORIGIN SPECIES	Human
ACCESSIONS	BC037955
VECTOR	pDONR223
NCBI SYMBOL	ADAT2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccaaagaagccctcgaaaatactgaagttctgttggcgtcttattggtctacaacaatgaagttgtaggggaaggggagaaatgaagtaaccaacacaaaatgctactcgacatgcagaaatgggtggccatcgatcaggtctctgattggtgtctcaaaagtggcaagagtcctct
gaagtattgaacacactgtgtgtatgtcactgtggagccgtgcattatgtgtgcagctgtctccgctgatgaaatcccgtggtgtgtatggctgtcagaatgaacgatttgggtgtgtgtggtctgttctaaatattgctctgtgacctacaaacactggagaccatttcagtgtatcc
ctggatcgggctgaggaagcagtggaatgttaagacctctacaacaagaaaatccaaatgcacaaaatcgaagttcggaaaaagggaatgtcagaaatctga