



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

CCR6-chemokine (C-C motif) receptor 6 Gene

Reference: PTXBC037960

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BN-1; C-C CKR-6; CC-CKR-6; CCR-6; CD196; CKR-L3; CKRL3; CMKBR6; DCR2; DRY6; GPR29; GPRCY4; STRL22; C-C chemokine receptor type 6; G protein-coupled receptor 29; LARC receptor; chemokine (C-C motif) receptor 6; chemokine (C-C) receptor 6; chemokine recepto		

For research use only.

GENE DESCRIPTION	chemokine (C-C motif) receptor 6
GENE ID	1235
ORIGIN SPECIES	Human
ACCESSIONS	BC037960
VECTOR	pDONR223
NCBI SYMBOL	CCR6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr