

TUG1-taurine upregulated 1 (non-protein coding) Gene

Reference: PTXBC037986

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CDCBM4; GCP-1; TUBG; TUBGCP1; tubulin gamma-1 chain; gamma-tubulin complex component 1; tubulin, gamma polypeptide; tubulin gamma 1		

For research use only.

GENE DESCRIPTION	taurine upregulated 1 (non-protein coding)
GENE ID	55000
ORIGIN SPECIES	Human
ACCESSIONS	BC037986
VECTOR	pDONR223
NCBI SYMBOL	TUG1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggactccagatttccagatttcaagacctggacctggaaccgaaagagctgtcacgatcgccgaggaacactggagtagattttttatttttaattttgggactgttgacctgtgtgagaaaagagacaacgactgagcaagcactaccaccagcactgttactggaattaga
agacctgagtttctgacagaccctcagtgcaactgaggatgctccatccaaagtgaattatgtcctgtgacctctgattgctgagtggttccactggacctctgactacctccctgtgtatttccatcagcctacagacctgtacctggatttttgcggagatgattctaccacttactactg
acgaagacacccattccagtggaacctgtga