



Reference: PTXBC037992

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GC-UNC45; GCUNC-45; GCUNC45; IRO039700; SMAP-1; SMAP1; UNC-45A; protein unc-45 homolog A; general cell UNC45; smooth muscle cell associated protein-1; unc-45 homolog A; unc-45 myosin chaperone A		

For research use only.

GENE DESCRIPTION	unc-45 homolog A (C. elegans)
GENE ID	55898
ORIGIN SPECIES	Human
ACCESSIONS	BC037992
VECTOR	pDONR223
NCBI SYMBOL	UNC45A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]