

MOBKL1A-MOB1, Mps One Binder kinase activator-like 1A (yeast) Gene

Reference: PTXBC038112

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MOBKL1A; MATS2; MOB4A; MOB kinase activator 1B; MOB1 Mps One Binder homolog B; MOB1, Mps One Binder kinase activator-like 1A; mob1 homolog 1A; mob1A; mps one binder kinase activator-like 1A		

For research use only.

GENE DESCRIPTION	MOB1, Mps One Binder kinase activator-like 1A (yeast)
GENE ID	92597
ORIGIN SPECIES	Human
ACCESSIONS	BC038112
VECTOR	pDONR223
NCBI SYMBOL	MOBKL1A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagctctgtttgtagtcgctcttctaaactttaaaccaagaagaacattccagaggggtctcaccagtagagctctttaaaccacgcagaagccacacttggcagtgccaacctcggatggctgcatgcttctgaaggggaagatctcaatgaatgggtgcagttaacactgtgatttctcaatcagatcaacatgctttatggaactatcacagactctgtacagaagagagttgtccagtgatgctcagctggccaaaatagatgatcattgggcagatggaacgaacataaagaacctaataagtgctctgcaccaagtatattgattactgatgactgggttcaggaccagttggatgatgagacgttattccatcaaaaattgggtcccttcccaagaattcatgctgtggcaaaaactatactcaaacgcctcttaggggttatgctcacattatcatcagcatttgaccctgtgatccagctcaggaggaagcacatcaaatcatcttcaagcactttatTTTTTgtccaggaattcaacttattgatagaagagaactgcaccactccaagaactgattgaaaactcactcaaaagacagataa