

CCL14-chemokine (C-C motif) ligand 14 Gene

Reference: PTXBC038289

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CC-1; CC-3; CKB1; HCC-1; HCC-1(1-74); HCC-1/HCC-3; HCC-3; MCIF; NCC-2; NCC2; SCYA14; SCYL2; SY14; C-C motif chemokine 14; chemokine (C-C motif) ligand 14; chemokine CC-1/CC-3; chemokine CC-3; hemofiltrate CC chemokine 1; new CC chemokine 2; small inducibl		

For research use only.

GENE DESCRIPTION	chemokine (C-C motif) ligand 14
GENE ID	6358
ORIGIN SPECIES	Human
ACCESSIONS	BC038289
VECTOR	pDONR223
NCBI SYMBOL	CCL14
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagatccgtggtgccattcccttctctcatcaccatgccttagggaccaagactgaatctctctcaggggaccttaccacccctcagagtgtcttcaactacactacactacaagatcccgcgtcagcgattatggattactatgagaccaacagccagtgctccaagc
ccggaattgtcttcatcaccaaaaggggccattccgtctgtaccaacccagtgacaagtgggtccaggactatatacaaggacatgaaggagaactga