



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

CCNL1-cyclin L1 Gene

Reference: PTXBC038394

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ANIA6A; BM-001; PRO1073; ania-6a; cyclin-L1; cyclin L ania-6a; cyclin L gamma; cyclin-L; cyclin L1		

For research use only.

GENE DESCRIPTION	cyclin L1
GENE ID	57018
ORIGIN SPECIES	Human
ACCESSIONS	BC038394
VECTOR	pDONR223
NCBI SYMBOL	CCNL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]