

C1orf105-chromosome 1 open reading frame 105 Gene

Reference: PTXBC038410

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C1orf105; chromosome 1 open reading frame 105		
For research use only.			

GENE DESCRIPTION	chromosome 1 open reading frame 105
GENE ID	92346
ORIGIN SPECIES	Human
ACCESSIONS	BC038410
VECTOR	pDONR223
NCBI SYMBOL	C1orf105
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaaaaagagaactaaaggcttctgttccaaaattgacaagattccttggttagtgaggccagccttgtaacaagccattagtgctcagcctccagaagatctctacacctctgcgactttctgacttcacgaagaagaatgaattgccaaattgtttcaagttccagatgttt
atctaaggccaggaggaaaccagtgtgactccatgctgctcagaacacacagctgtgctccacatgcaagaatgaaatggtaaccaagaacaatgaaatccagatgatccaaagcatcctttgagaattgtatgagttatagaatgagcttcatcaacccaattccagactaca
cctgagcctttcatgatgacatcccaacagaaaacttactacagactgccattctgggcccaggacagctgtctccacggattactgacagaggcctacaaaactctaaaagagagacaacgttcttcttgcagaaaggaaccaataggcaagacaacgagcagtgga