

TSGA13-testis specific, 13 Gene

Reference: PTXBC038415

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	testis-specific gene 13 protein; testis-specific A13; testis specific 13		
For research use only.			

GENE DESCRIPTION	testis specific, 13
GENE ID	114960
ORIGIN SPECIES	Human
ACCESSIONS	BC038415
VECTOR	pDONR223
NCBI SYMBOL	TSGA13
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagccaaaagagacaaacaaagtttcaaaatggcaaatcaagacttcagaaaatagctcagctaaacgtgagaaaggaaatggtgtcaatagcaaaagattctgatgcagtcgggcaatcaaaattgtctagagaacctcggcattacagtcctccaaattggccagctact
ataagcctttgaaggccactgcctgcagaaattcctggctcaaaacaggaaaaacaccagcttcatgtttaaagtaacacagtagtaccaggataagacgttactgattatgaccaacaacccacctcctgctcaatccccagcaagacaaggagagtcacataatattttccaaagg
agttactgctcaaggatgaaagtcacatcagcacaacccactgagaacctctgctgccccgcagctcagaaaaaaagttaagatctaagctgaaaccaatcttcccttggatgactgctggatgatccacatccaagcgagaacaatggttaggtttccactgacaatgattca
agagcgaagggaagtattcaaaagtctacgtcttggagacacagaagaaatgtaccctcagctcacttggctcagtcctcagtgagagagatatgaggaaagatgcttcaagaagtcagcagtgaaaggccaattccaagtattcgggaaccactgacactgcctccttggga
agacatgccaccagaaccgcggcgaggagagcgcttcgcaacgggagggccccgcagtggtattatcaaaaggccacgggtcctcgggtga