

SNF8-SNF8, ESCRT-II complex subunit, homolog (S. cerevisiae) Gene

Reference: PTXBC038830

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SNF8, ESCRT-II complex subunit; SNF8, ESCRT-II complex subunit, homolog; vacuolar-sorting protein SNF8; Dot3; EAP30; VPS22; EAP30 subunit of ELL complex; ELL-associated protein of 30 kDa; ESCRT-II complex subunit VPS22		

For research use only.

GENE DESCRIPTION	SNF8, ESCRT-II complex subunit, homolog (S. cerevisiae)
GENE ID	11267
ORIGIN SPECIES	Human
ACCESSIONS	BC038830
VECTOR	pDONR223
NCBI SYMBOL	SNF8
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcaccgcgcgggtggagctggcgcctcgcgaagaagaacttcagagggccaagtataaggagcgaggagcggctctggctgaggaccagctagccagatgtcaaagcagttggacatgttcaagaccaacctggaggaaattgccagcaaacacaagcaggagatcc
ggaagaatccctgagttccgtgtgcagttccaggacatgtgtgcaaccattggcgtgcatccgtggtctctggaaggaattttgtctgagatgctggcgtgggggacttctattacgaactaggtgtccaaattatcgaagtgtgcctggcgtgaagcatcggaatggaggtctgataa
ctttggaggaaactacatcaacaggtgtgaaagggaagggcgaagtgcggcagatgtcagtcagatgacatgacagagccatcaagaactaaagcacttggcactggctcggcatcatccctgtggcggcacttacctcattcagctgttcagctgagctcaatagatcaca
ccgtgtgtcagctggcagagaatggctacgtgactgtcagtgagatcaagccagctcttaaatgggagaccgagcgagcggcgaagtgtctggaacacctgctgaaggagggttggcgtggctggacttacaggccccaggggaggccactactggctgccagctcttcac
tgacctctactccaggagattacagctgaggaggcagagaagccctccctga