

PRR16-proline rich 16 Gene

Reference: PTXBC038838

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DSC54; LARGEN; protein Largen; mesenchymal stem cell protein DSC54; proline-rich protein 16; proline rich 16		

For research use only.

GENE DESCRIPTION	proline rich 16
GENE ID	51334
ORIGIN SPECIES	Human
ACCESSIONS	BC038838
VECTOR	pDONR223
NCBI SYMBOL	PRR16
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgactgacagctccaaacggacacgctgaatagtagctcaagtgacacaacagcctccagcctagagaagatcaaatgacagctaatgcaccgttataaacccccagcacacacatctgtatctcaggtcctgagaagccaaacctccaccctctccaaggtgacacctgtgaaagtfgaagaccccaaaagggtgtcccaactccaatctgtfaaaaaccaatggcacctcttacgaaatggagcttaccaggtggacctaacaaaattccaaatggagatctgtgcatacccaacagtaactggacaaggctccagcttctgatcctagacctgaaaaagacagatgtccccagcaggcctcgagaacgagttcggtttaataaaaaagtacagctacatggctattgtctgactgtgataccgggtataacataaaaaacagggaggtccacttacacagtgaaactgtccaccaccgggaagattcctaccaaggccctccctccctctacacccccctctccctcttcccactagaaaatggggaaatgggaataagccacagtaacagcttccccctatcagacctgcaactgtcctctcccactgcacaaaaccacagaagacgatctgaggagtaaccactacaaccgtgtga