



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

NAALAD2-N-acetylated alpha-linked acidic dipeptidase 2 Gene

Reference: PTXBC038840

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GCPIII; GPCIII; NAADALASE2; NAALADASE2; N-acetylated-alpha-linked acidic dipeptidase 2; N-acetylated alpha-linked acidic dipeptidase II; NAALADase II; glutamate carboxypeptidase III; N-acetylated alpha-linked acidic dipeptidase 2		

For research use only.

GENE DESCRIPTION	N-acetylated alpha-linked acidic dipeptidase 2
GENE ID	10003
ORIGIN SPECIES	Human
ACCESSIONS	BC038840
VECTOR	pDONR223
NCBI SYMBOL	NAALAD2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr