

C1orf91-chromosome 1 open reading frame 91 Gene

Reference: PTXBC038842

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C1orf91; AASL548; PRO1105; RP4-622L5; dJ622L5.7; transmembrane protein 234		
For research use only.			

GENE DESCRIPTION	chromosome 1 open reading frame 91
GENE ID	56063
ORIGIN SPECIES	Human
ACCESSIONS	BC038842
VECTOR	pDONR223
NCBI SYMBOL	C1orf91
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggcgtctctggggcaggtgtggctctggctgtggcgcgtctgtggcgcagccgctgtgaagcggcctccgcccgcctgcagcgggttcagagcgcagcctggggccagcagttgctacaggagatgaagaccctctcttgaatactgagtacctgatgcc
ctttctctcaaccagtgtggatccctctctattacctcacctggcatcgacagatctgacctggctgtgcccactctgtaactctctggctatcatcttcacactgattgttgggaaggccctggagaagatattggtggaaaacgagcagttgctgcatggtgctcaccgtgataggaaattc
actctgcatcacaagctcagttccatggactgcagaactccagctgcatggaaggccagctgcagactttgagccagaaatgcaaacggaggccctctgggactcagtcagagcgtttggctga