

LY86-lymphocyte antigen 86 Gene

Reference: PTXBC038846

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MD-1; MD1; MMD-1; dJ80N2.1; lymphocyte antigen 86; MD-1, RP105-associated; ly-86; protein MD-1		

For research use only.

GENE DESCRIPTION	lymphocyte antigen 86
GENE ID	9450
ORIGIN SPECIES	Human
ACCESSIONS	BC038846
VECTOR	pDONR223
NCBI SYMBOL	LY86
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagggtttcacagccactctctctctggactctgattttccagctgcagtgaggcggcgggtggaaagcctggccacacacgtggtctgtagcgacagcggcgttgaagtgcctaccagagttgcgatccattacaagattttggcttttctgtgaaaagtgttccaagcaatta
aaatcaaatatcaacattagatttgaattattctgagagagacatcaaaagcgttttttgacctagctctcatgtctcaaggctcatctgttttgaatttctctatcccatctgtgaggcggcctctgcccaagttttcttctgtggaagaaggaaaggagagcagatttactatgctggcctgt
caataatcctgaatttactatctcaggaggagaataccaggttttctggaactgtacactgaaaaacgggtccaccgtggcctgtgccaatgctactatcatgtgtcctga