

TMEM126B-transmembrane protein 126B Gene

Reference: PTXBC038933

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	complex I assembly factor TMEM126B, mitochondrial; HT007; transmembrane protein 126B		

For research use only.

GENE DESCRIPTION	transmembrane protein 126B
GENE ID	55863
ORIGIN SPECIES	Human
ACCESSIONS	BC038933
VECTOR	pDONR223
NCBI SYMBOL	TMEM126B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagcatctatgcatggcagccagtcctctctagaagatgcaaaactcagaagacaaatggcagaaatcatagaaaaaatttgactatcttagaaaaaagaatgacacaaaatatatacaaatggcgacattggaacaacagctggttctctggaatattctcaaactcctgtt
cagacgctgctcaaggttaaacatgagctttgaagacatgcatcattggctacacttccattttgtctactgtgttactgacaaagctttttgtaattgatgctttgtatcagataataagaagcaagaaaaactgtgtttcagaagctcactgatggcatagtttgtgtgttttctatccagttctt
tggcttttactaaaaatggagcctggcaaccaagatcataccgttccactgccacaaaaggagggttttaaccattggatgacgctttgtcaacacaaatgaaatgaatggcgattcctctagcttctcagattatgtttggaatattaatgggtctataccattatgcagattttgaagagac
acttgagaaaactatacatgaagagtaa