

ANKRD19-ankyrin repeat domain 19 Gene

Reference: PTXBC038951

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	ankyrin repeat domain 19
GENE ID	138649
ORIGIN SPECIES	Human
ACCESSIONS	BC038951
VECTOR	pDONR223
NCBI SYMBOL	ANKRD19
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggaagctctcagcttgggagagcctgggagggcgtcctggactccatggaccaagagtaacgtgctgggggtaccacatcgggactgggaactgcggaagatccacagggtgcatcaaggcgacgcccgggaggtggagcactgcctgacgcgaggttcc
gggacttggagcccgcgacagaaaagacaggactgttctacatttgacctgtgcccattggccgtgtggaagtgtgacccctctgtgagcagaaatgccagatcaacatatatagactaaacaggacacctttaatgaaggctgtacactgccaggagaggttctgtccattattc
tcctggaacatggcgccaatccaaacatttaaggatatctacagcaacactgctcctattatgctgtgtataataaggggacttctcagagaaaaactgcttcccaccatgcaaatattgaagcactaaacaggagggaacacatccactttgttctataaattctaggagacagcaaa
ttgtggaatttctgttgaagaaccaggcaaatftacatgccattgacaatttcagaaagacagccctcatgcttgcagtacagcataactcatcaagtatcgtcagcctctcctcaacaaaataataatcttttctcaagacctgtttggcagactccgaggattatgctgtttgtataattca
gaagcattcaacaacaattttggaacataaaaaataataacttaaaagtcattcttga