

## EHF-ets homologous factor Gene

Reference: PTXBC038995

|                |                                                                                                                                                                       |                     |                |
|----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                   | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                  | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | ESE3; ESE3B; ESEJ; ETS homologous factor; ESE3 transcription factor; ETS domain-containing transcription factor; epithelium-specific Ets transcription factor 3; hEHF |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | ets homologous factor                                                                             |
| <b>GENE ID</b>          | 26298                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC038995                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | EHF                                                                                               |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atgattctggaaggaggtgggtgaatgaatcaacccggcaacaacctctccaccagccgagcctggacagacagctactccacgtgcaatgtttccagtggtttttggaggccagtgatgaaattcatcctcagtagtgaccaagtaccaggtgtgggagtggtccagca  
cctcctggacaccaaccagctggatgccaattgtatccctttccaagagttcgacatcaacggcgagcacctctgcagcatgagtttcaggagttcaccggcgccaggagacggcgggcagctcctctacagcaacttgagcatctgaagtggaaaggccagtgagtagtgacct  
gttcaggtccacacacaaatgtcattgtcaagactgaacaaactgagccttccatcatgaacacctggaagagcagagaactattatgacaccaactatggtgacacagtagatttgggacagcaaaactttctgccgggctcagatcctatgacaaccaccagtcacctctgttgacag  
agtcacctgatatgaaaaaggagcgaagaccctcgccaagtggcacacacaaaaagcacaaccggagaggggactcacttatgggaattcatccgcgacatcctctgaaccagacagaagaccaggatgattataaaatgggaagaccgagctgagggcgctctcaggttctgaaatcag  
aggcagtggtcagctatgggtaaaaagaagaacacagcagcatgacatgaaaagctcagccgagctatgagatattactacaaaagagaattctggagcgtgtggatggacgaagactggtatataaattgggaagaatgccgaggtggagagaaaatgaaactga