



Reference: PTXBC039015

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LRH-1; P2X5; P2X5R; P2X purinoceptor 5; ATP receptor subunit; ionotropic ATP receptor P2X5; lymphoid-restricted histocompatibility antigen-1; purinergic receptor P2X, ligand gated ion channel, 5; purinergic receptor P2X 5		

For research use only.

GENE DESCRIPTION	purinergic receptor P2X, ligand-gated ion channel, 5
GENE ID	5026
ORIGIN SPECIES	Human
ACCESSIONS	BC039015
VECTOR	pDONR223
NCBI SYMBOL	P2RX5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]