



Reference: PTXBC039026

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	fas-activated serine/threonine kinase; FAST kinase; Fas activated serine/threonine kinase		

For research use only.

GENE DESCRIPTION	Fas-activated serine/threonine kinase
GENE ID	10922
ORIGIN SPECIES	Human
ACCESSIONS	BC039026
VECTOR	pDONR223
NCBI SYMBOL	FASTK
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]