



Reference: PTXBC039071

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein-cysteine N-palmitoyltransferase HHAT; MART2; SKI1; Skn; melanoma antigen recognized by T-cells 2; skinny hedgehog protein 1; hedgehog acyltransferase		

For research use only.

GENE DESCRIPTION	hedgehog acyltransferase
GENE ID	55733
ORIGIN SPECIES	Human
ACCESSIONS	BC039071
VECTOR	pDONR223
NCBI SYMBOL	HHAT
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1