



Reference: PTXBC039148

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	K6L; KRT6L; keratin, type II cytoskeletal 79; CK-79; K79; cytokeratin-79; keratin 6L; keratin 79, type II; type-II keratin Kb38; keratin 79		

For research use only.

GENE DESCRIPTION	keratin 79
GENE ID	338785
ORIGIN SPECIES	Human
ACCESSIONS	BC039148
VECTOR	pDONR223
NCBI SYMBOL	KRT79
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr