

WFDC5-WAP four-disulfide core domain 5 Gene

Reference: PTXBC039173

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PRG5; WAP1; dJ211D12.5; WAP four-disulfide core domain protein 5; p53-responsive gene 5 protein; protease inhibitor WAP1; WAP four-disulfide core domain 5		

For research use only.

GENE DESCRIPTION	WAP four-disulfide core domain 5
GENE ID	149708
ORIGIN SPECIES	Human
ACCESSIONS	BC039173
VECTOR	pDONR223
NCBI SYMBOL	WFDC5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggaccagagccttctcctctgggggcccctctggtgtgtgggagtcagctgcctgtgtctttggcaggaagaaggagagaaatcggggggctgcccgccagatgatggccctgcctctctatcgggtgcctgaccagtgctgtgaagacagccagtgcccttgaccagga
agtgtgtctacagagctgtgtccgagtggtgtcccgagggtctgtgtgaagctgggcagctgccagaggaccaactgcgtgcctcagcccatgaaccacctgtgtcacaaggactcagactgctcgggcaaaaagcgatgctccacagcgcctgcggcggtgattgccggga
tcctgccagaggctaa