



Reference: PTXBC039174

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CMD1W; CMH15; HEL114; MVCL; epididymis luminal protein 114; meta-vinculin; metavinculin		

For research use only.

GENE DESCRIPTION	vinculin
GENE ID	7414
ORIGIN SPECIES	Human
ACCESSIONS	BC039174
VECTOR	pDONR223
NCBI SYMBOL	VCL
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr