

TMEM47-transmembrane protein 47 Gene

Reference: PTXBC039242

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BCMP1; VAB-9; transmembrane protein 47; brain cell membrane protein 1; transmembrane 4 superfamily member 10		

For research use only.

GENE DESCRIPTION	transmembrane protein 47
GENE ID	83604
ORIGIN SPECIES	Human
ACCESSIONS	BC039242
VECTOR	pDONR223
NCBI SYMBOL	TMEM47
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcttcggcggcagcggcatggagggtgcgcgtgtcgtgtgaccccttgaagctgtcggcgtgtgtgcatcttcctggcgtgtgtgacctggggcggtgctgagcccgctgggtcacagctgaccaccagtactacgtgtgtgtggagtcctgccgcaa
acccgccagcttggacatctggcactgcgagtcacgcagcagcgaattggcagattgctactctgcttactcctggcggcgctgccatcatctcattgcattcctgtgggttgattctatctgcgtgggacatcgaaggcgtttctatagacctgttgcggcatgcttttgcagcagt
tgttttacaggtttgcagcctgggtcctttacccaatcaagttcattgaaactgtgagctgaaaatttaccatgagttcaactgggggttatggcctggcctgggggtgcaactatatttctgtttgggggtgccatcctttatgacctgaaccctaagaactatgaagactactactag