

## CPA5-carboxypeptidase A5 Gene

Reference: PTXBC039362

|                               |                                                      |                     |                |
|-------------------------------|------------------------------------------------------|---------------------|----------------|
| Product type:                 | DNA                                                  | Storage Condition:  | -20°C          |
| Size/Quantity:                | 2 ug                                                 | Delivery Condition: | +4°C, Blue Ice |
| Synonyms                      | carboxypeptidase A5; testicular tissue protein Li 32 |                     |                |
| <b>For research use only.</b> |                                                      |                     |                |

---

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | carboxypeptidase A5                                                                              |
| <b>GENE ID</b>          | 93979                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC039362                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | CPA5                                                                                             |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgcaggcactctggaggcgggacgcgcctgggccatccccctgggacagcgggacacactcctgttctcagctttatcctggcagcagcttgggccaatgaatttcacaggggaccaggtcttcgagctctggccaagatgagaagcagcttcacttctcgggatctggag  
ggcctgaaacccagaaaggtgacttctggcgtggccagccaggcccgctccctgtggatagagagttccttctccgaactgaaagacatcaaaagcttatctggagctcatggacttctacagcatcatgataaaggacatccaggtgctgctgaggggaagacagccat  
ggcgaatcccgcggctggagcgcagcaccaacagcttcagttactcatcacacacccctggaggagatatagctggattgacaacttgaatggagcattccgataattgtctcaaaatcagattggcaacagcttggaaaccagtcattcttgcctgaagttagcactggag  
gttctcggcaccagccatctggatcgacactggaaftcactccccggagtgatcacccatgccacggcatctggactgccaataagattgctagtgattggcaaaagaccgtgtcctgacagacatactgaatgccatggacatcttcatagagctcgtcacaacccatgatgggttg  
ctttaccacagcatgaaccgttatggcggaagaacaagtcacagactggaaatctctgcatcggtggatctcaacaggaactggaatcggttttggaggaaatggttcaacagcaacccctgctcagaacttatcacggccctccctcagtcggagcggaggtggct  
gccatagtgaactcatcacagcccatggcaactcaaggtctgactcctccacagctactctcagatgcttatgaccttacggccgattgctggagcccgtttcaaatcagaggaggtgtacgatcttccaagatcggtggagccctgtataaggtccatgggatcagatcat  
ttttggcagcatcagcaccacctctatgtggccagtggtgacacgtcgaactgggctatgacagtggtcatcaagtacgccttcagctttagctccgggacactgggcagtatggcttctgctgcccggccacacagatcatccccacggccaggaagcgtggtggcgttcggacc  
atcatggagcacacctgaatcacccctactag