

LOC730441-similar to trypsin X5 Gene

Reference: PTXBC039387

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	similar to trypsin X5
GENE ID	730441
ORIGIN SPECIES	Human
ACCESSIONS	BC039387
VECTOR	pDONR223
NCBI SYMBOL	LOC730441
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtctatctcagtcctccagaacctgcgtgggtctctcattcaccttgactgggtattgacagctgccactgcccttacctgttgaaattcgactgggagttctcaacctagcatcacaacaagaacaacagatatga