

C6orf66-chromosome 6 open reading frame 66 Gene

Reference: PTXBC039464

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C6orf66; HSPC125; My013; bA22L21.1; NADH dehydrogenase [ubiquinone] 1 alpha subcomplex assembly factor 4; NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, assembly factor 4; NADH dehydrogenase (ubiquinone) complex I, assembly factor 4; hormone-regulate		

For research use only.

GENE DESCRIPTION	chromosome 6 open reading frame 66
GENE ID	29078
ORIGIN SPECIES	Human
ACCESSIONS	BC039464
VECTOR	pDONR223
NCBI SYMBOL	C6orf66
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggagcactagtgtattcgggtatcaggaattcaacctagagaaccgagcggaacgggaaatcagcaagatgaagccctctgtcgtccagacacccctctaccaacagccctctgcgagagcagattagtctctatccagaagttaaaggagagattgctcgtaaagatgaaag
ctgctgtcgtttctaaaagatgtgtatgttattccaaagatcctgtgtcttccttgcaggtaaaagctgctgaacatgtcaagagccgaaggaaatcagattgccgaagaccatcatttgatatgataaataatgaagcattcccaaggcaaaattccattgtagaagcattgacactctca
ataatcataagcctttccagaaacctggactgctgagaaaataatgcaagaataccagctagaacagaaagatgtgaattctctttaaattttgtacttttgaagtcgaaatctccctctgaagacaagaagcaatcacgatacaaatga