



Reference: PTXBC039574

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MIG1; SKD1; SKD1B; VPS4-2; vacuolar protein sorting-associated protein 4B; cell migration-inducing 1; cell migration-inducing gene 1 protein; suppressor of K(+) transport growth defect 1; suppressor of K+ transport defect 1; vacuolar protein sorting 4B; v		

For research use only.

GENE DESCRIPTION	vacuolar protein sorting 4 homolog B (S. cerevisiae)
GENE ID	9525
ORIGIN SPECIES	Human
ACCESSIONS	BC039574
VECTOR	pDONR223
NCBI SYMBOL	VPS4B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]