

## FAM43A-family with sequence similarity 43, member A Gene

Reference: PTXBC039689

|                               |                                                             |                     |                |
|-------------------------------|-------------------------------------------------------------|---------------------|----------------|
| Product type:                 | DNA                                                         | Storage Condition:  | -20°C          |
| Size/Quantity:                | 2 ug                                                        | Delivery Condition: | +4°C, Blue Ice |
| Synonyms                      | protein FAM43A; family with sequence similarity 43 member A |                     |                |
| <b>For research use only.</b> |                                                             |                     |                |

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | family with sequence similarity 43, member A                                                     |
| <b>GENE ID</b>          | 131583                                                                                           |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC039689                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | FAM43A                                                                                           |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgctgccgtggaagaagcacaagttcagctgctggccgagggccgcccggcaggcgtccaagcccaaggctacgtgtgagcctgcactactcggcgtcagctcgtggcgcggcggtccccgaaggcgcgttagccgggtgggcagcatgttcgctccaagcgca  
agaagctgcacatcactagcaggaccacacttacaccgtgctctacctgggcaatgccaccacatccaggcgcggcgacggctgcaccgacctgtgtgggcaagatctggagcaagagcaggcggccgtcaggccaccaagatgaagctgacgtgagtgccgaggg  
tatccgcatggtgcacgcgagggagcgcgctgcccggccgcccggccacctctacctgctgcaccgctgacactactgctgcccgcgacgcggctgcccgaaggtcttcgctgggtgtaccgacacgagctgaagcacaaggccgtgatgctgctgccacgctgctggtgtc  
caagccgaaaaaggcagggccatggccctgctgctctaccagacgtcggccaacgcgtggcggaattaaacgcctcaagcggcgagcgcgctcaccagcagcaggagctgggtggcgacacacacatcccgtagtgccgctgcgcaagctgctctacacggacc  
ctgctgtataaacgccgggtggagcgcagcgcgagcgcgcccagctggctccatcacaggagacctgctcgggaacagctggagcaggagctgcaggaggaagaggaagagcaacccgagggctgcccggaggaaggaggaacgtgcggcagagggagatcc  
agcagaggaggagccgagggcgcagcgtgctgtagtggtgcacatgcactttgagtgccgggactgtgtgatactctggagaatggccgtggggagggcctaggagggcggcggtccctggcccgggggccggccgcccctctgctgctggcagcgcctccgacat  
gaagctgagctgtcgaacttattagcgacctggggcagctcagcttggcaacgacgtgcacccctgcaggccgacttgcgggtgacgcgcctgctgtcggcgacagcacgggcagcgagagctccatcgagggcgggggccctgacgccacctccgaccgcccggga  
ctgtccccccaggccgacggcgccagtgcagacgagccccactcgggctga