

LOC339047-hypothetical protein LOC339047 Gene

Reference: PTXBC039707

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical protein LOC339047
GENE ID	339047
ORIGIN SPECIES	Human
ACCESSIONS	BC039707
VECTOR	pDONR223
NCBI SYMBOL	LOC339047
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttttctgcttaggatgaatggctgagcggaggctgtaaaacctggcactctgcttgggttatcaatactctgctgacctgctatcgtgggactgacttgggtggaagtccttgggttacttacttactgtgtttctgagaagtataaattgccatccctctgcacaagttacctttgtg
tgtcttctgaagactatcttccgtctcaaatggacatgatggatccacggatgtacagcagagagccaggggtccaactgccgtagacaggaaggaaattaaattgtcctggagacatcttactttatggagacaggttgaaaccaaagttcgagctaaatccgtaagatgaaggt
gacaacaaaagtcaaccgtcatgacaaaatcaatggaaagagggaagaccgcaagaacatctgaggaactaagcatgaaagaacgtgagcagggagaaaaaggagggcaggtgtcagaggcagaggaaaatgggaattggatatgaaagaatacacactacatggaatg
ttcaacgtgcgcaagcgttgcggcgggcagaggactactacagatgcaaaatcacccttctgcaagaagccttcttgaaccgggtgcacttacaattgtgtattactctggagtgaaaagaataattacagccatgcctaa