

ANKRD10-ankyrin repeat domain 10 Gene

Reference: PTXBC039715

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ankyrin repeat domain-containing protein 10; ankyrin repeat domain 10		
For research use only.			

GENE DESCRIPTION	ankyrin repeat domain 10
GENE ID	55608
ORIGIN SPECIES	Human
ACCESSIONS	BC039715
VECTOR	pDONR223
NCBI SYMBOL	ANKRD10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcggcgccgggagcggcgccggcgtagggcggtctccagcggaggagctgctctgctcgtttcccgctgcacgcgcctgcgcgcgacggcgacctggccacgctctgctcgtcgcagcagacacccacgcccacctggcctctgaggactccttctatggctgg
acgcccgtgcactggcgccgcatctcggcaagttggagtgcttagtgagtgagagcgggagccacactcaacgtctccacacacggtagcgcgcagacgccacattgcagcctttgggggacatcctcagtgccctggtctggtgattcaagcaggagccaacattaa
caaacgggattgtgagggtgaaactccattcacaaggcagctcgtctgggagcctagaatgcatcagtgcccttggcggaatggggtcacgtcgataatcctagaagggggagcattccatag