

PPHLN1-periphilin 1 Gene

Reference: PTXBC039832

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HSPC206; HSPC232; periphilin-1; CDC7 expression repressor; gastric cancer antigen Ga50; periphilin 1		

For research use only.

GENE DESCRIPTION	periphilin 1
GENE ID	51535
ORIGIN SPECIES	Human
ACCESSIONS	BC039832
VECTOR	pDONR223
NCBI SYMBOL	PPHLN1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcttacagaagagacgaaatgtggtctgagggacgatgatgaatgaagaattccgagagaacgagcacctctcgaagtcacccagtgatgaatctggttatagatggacaagagacgatctctgcaagcaggcaacctgaatacaggacatgagagatggctttaagaaga
aaaagtgttctactctccattatgcgagagagcgggtctctataaaagggaacaatactttttcagagaatcacctgttgccgaaaggattctccacacagcagatctggttccagtgtagcagaagctactctccagaaggagcaaatcactcttccatcagtcacatagaa
agtccgtgctgctggtgctcctacaacggcagaatgaagaaatcctgaagagataaagagagggcctgtccagctcttgaacacatcaagagatacttcaccccaagtggttcagcagtttctcatcaaaaggtgttagacaacccagtaggctaactgaaaaggaaactgtctgagg
ctgcaagcaagtggtgctgctgaaaagctagagaaatcagatgaaagtaacttgcctgaaatttctgagtatgagggggatccacagcaccattgttactgaccagccagaggaaacctgagtcacaacaacacatgggtagaattatttgaagatagtcagtaaccactgcctctaa
gcaatagcatcaaaaaccaaagagattgaacagggtgagagctag