



Reference: PTXBC039859

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane and coiled-coil domains protein 1; transmembrane and coiled-coil domain family 1		

For research use only.

GENE DESCRIPTION	transmembrane and coiled-coil domain family 1
GENE ID	23023
ORIGIN SPECIES	Human
ACCESSIONS	BC039859
VECTOR	pDONR223
NCBI SYMBOL	TMCC1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr