

ODF3L1-outer dense fiber of sperm tails 3-like 1 Gene

Reference: PTXBC039862

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	outer dense fiber protein 3-like protein 1; outer dense fiber of sperm tails 3 like 1		
For research use only.			

GENE DESCRIPTION	outer dense fiber of sperm tails 3-like 1
GENE ID	161753
ORIGIN SPECIES	Human
ACCESSIONS	BC039862
VECTOR	pDONR223
NCBI SYMBOL	ODF3L1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaactgccaaagggaccaggagctctgtgtactttgcacagcaccagaaaaggagccattgccctcaaggcaggaggtcaagcagacccctgtcatcatggccaagatcaaaggtccggggcccgccaagtacctccggccatctgcacgggtacatagatcatgacatct
ccatgttcaaggcaccagcttatacctgcatagccggcactcagagaagcggatggtgtgcccacagcagccctgggcttgcctctcttgatcccaataactcgggttggaatgtccagctgcccgcaggtcccatggaggagcgcctccaacctgcgctgaacccccacct
cgcatcctgccagtactctttgagaagatccaccaccgggggaacgcagggtctcccgatcacgtttggctaccggcgcccatcacagagtgtggacctcaacccgggtcccaaccagtaccagatgccactcttctggggcccaacacctgtcagccgagctgctccctgct
acagtctggcctccagggacaagaactggttctacaaggagatgtgacagggagccctggacctaccaggtacggccgacctgagccatccatctatcagaaccgcagccctacttacagatggccaagcgttcgctaccctctggacctcacgcccacggcctggcccggtc
ccagaggtccagcaggtcactgtgcacaagcccacatccctgtttaccatgggcatcaagcactcactccacctgtgccactggtcatcgacattcgtgactga