



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

NRM-nurim (nuclear envelope membrane protein) Gene

Reference: PTXBC039865

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NRM29; nuclear rim protein; nurim (nuclear envelope membrane protein)		
For research use only.			

GENE DESCRIPTION	nurim (nuclear envelope membrane protein)
GENE ID	11270
ORIGIN SPECIES	Human
ACCESSIONS	BC039865
VECTOR	pDONR223
NCBI SYMBOL	NRM
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]